

IHE Integration Statement			
Vendor	Product Name	Version	Date
Clinisys	Clinisys GLIMS	9.8	4/12/19

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Application	None
Consistent Time	Time Client	None
Laboratory Code Set Distribution	Code Set Master	None
Laboratory Testing Workflow	Department System Scheduler/Order Filler	Report Fac-Simile For Order Group
Laboratory Testing Workflow	Department System Scheduler/Order Filler	None
Patient Administration Management	Patient Demographics Consumer	Merge
Patient Administration Management	Patient Demographics Consumer	None
Patient Administration Management	Patient Encounter Consumer	Maintain Demographics
Patient Administration Management	Patient Encounter Consumer	None
Sharing Lab Report	Content Creator	None
Cross-Enterprise Document Sharing	Document Source	None

Internet address for vendor's IHE information :

<https://www.clinisys.com/in/en/products/clinisys-glims-en/>

Links to Standards Conformance Statements for the Implementation

<b>HL7</b>	<a href="https://www.hl7.org/">https://www.hl7.org/</a>
<b>DICOM</b>	

Link to general information on IHE : [www.ihe.net](http://www.ihe.net)